

TCM5 Connected devices

Updated: August 2023



	Manufacturer	Device	Connectivity available	Manufacturer SW version	From TCM5 SW version	Comments	
Patient Monitoring Systems	Dräger	Infinity Delta	✓	VF10.1-W	1.3.0	Only with Sensor 84	
		Infinity Delta XL	✓	VF10.1-W	1.3.0	Only with Sensor 84	
		Infinity Kappa	✓	VF10.1-W	1.3.0	Only with Sensor 84	
		Infinity Omega	✓	VF10.1-W	1.3.0	Only with Sensor 84	
	Mindray	Beneview T5 / T6 / T8 / T9	✓		1.2.0		
	Philips	IntelliView MP40-90	✓			1.1.8	Using Philips EC10 Interface Module or IntelliBridge EC10 integral interface board with Open Interface Driver (ED101) and EC5 ID Module #102
		IntelliView MX400-800	✓			1.1.8	Using Philips EC10 Interface Module or IntelliBridge EC10 integral interface board with Open Interface Driver (ED101) and EC5 ID Module #102
	GE Healthcare	Carescape B450	✓	3.1.353 6ICU and 2 6ORP		1.3.0	
		Carescape B650	✓	3.1.353 6ICU and 2 6ORP		1.3.0	
		Carescape B850	✓	3.1.353 6ICU and 2 6ORP		1.3.0	
		Dash 5000	✓	7.3 MULTI 26MAR15		1.3.0	
	Nihon Kohden	BSM series	✓	V.01-03		1.2.0	
		CSM series	✓	Ver.00-08 12		1.4	
	Fukuda Denshi	DYNASCOPE 8000 Series	✓	V96-02		1.4	
Polysomnographs	Cidelec	CID analog	✓			1.2.0	
	NOX Medical	A1	✓	Noxturnal software version 5.1.0.19071		1.4.0	Connection available through the NOX C1 wireless access point
	Compumedics	E-Series	✓			1.3.0	Compatible with ProFusion Sleep software, version 4.2 or later
		Grael	✓			1.3.0	Compatible with ProFusion Sleep software, version 4.2 or later
		Siesta	✓			1.3.0	Compatible with ProFusion Sleep software, version 4.2 or later
Somt�		✓			1.3.0	Compatible with ProFusion Sleep software, version 4.2 or later	



	Manufacturer	Device	Connectivity available	Manufacturer SW version	From TCM5 SW version	Comments
Polysomnographs	Natus	Embletta Gold	✓		1.2.0	
		Embla S7000 / N7000	✓		1.2.0	
		Embla S4000	✓		1.2.0	
		Embla S4500	✓		1.2.0	
		Mdrive	✓		1.2.0	
	Philips Respironics	Alice analog	✓		1.2.0	
	Löwenstein	Miniscreen PRO	✓		1.3.0	
Ventilators	Breas	Nippy4+	✓		1.3.0	
		Vivo45	✓		1.3.0	
		Vivo45LS	✓		1.3.0	
	eVent Technologies	eVolution	✓		1.4.0	
EMR	Capsule Technologies	Medical Device Information Platform	✓		1.3.0	
CPET	Cosmed	Quark RMR	✓		1.4.0	
Data analysis SW	Stowood	Visi-Download	✓		1.4.0	

Legend: ✓ Connection validated